

## Datasheet

### Nbn polyclonal antibody

**Catalog Number:** PAB15515

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of Nbn.

**Immunogen:** A synthetic peptide corresponding to amino acids 150-200 of mouse Nbn.

**Host:** Rabbit

**Reactivity:** Human, Mouse

**Applications:** IHC-P, WB-Ce

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This antibody is useful for Western blot, where a band is seen at approx. 97 KDa.

**Form:** Liquid

**Recommend Usage:** Western Blot (0.5 ug/mL)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:100)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (30% glycerol, 0.1% sodium azide)

**Storage Instruction:** Store at 4°C for short term. For long term storage store at -20°C.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 27354

**Gene Symbol:** Nbn

**Gene Alias:** Nbs1