

Datasheet

EIF2B5 mouse monoclonal antibody (hybridoma)

Gene Symbol: EIF2B5

Gene Alias: CACH, CLE, EIF-2B, EIF2Bepsilon, LVWM

Catalog Number: H00008893-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant EIF2B5.

Immunogen: EIF2B5 (AAH13590.1, 1 a.a. ~ 721 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MAAPVVAPPGVVSRANKRSGAGPGGSGGGGARGA
EEEEPPPLQAVLVADSFDRRFFPISKDQPRVLLPLANV
ALIDYTLEFLTATGVQETVFCCWKAQIKEHLLKSKW
CRPTSLNVVRIITSELYRSLGDVLRDVKALVRSDFLL
VYGDVISNINITRALEEHRLRRKLEKNVSVMTMIFKESS
PSHPTRCHEDNVVAVDSTTNRLVLFQKTQGLRRFAF
PLSLFQGSDDGVEVRYDLLDCHISICSPQVAQLFTDNF
DYQTRDDFVRGLLVNEEILGNQIHMHVTAKEYGARVS
NLHMYSAVCADVIRRWVYPLTPEANFTDSTTQSCTHS
RHNIYRGPEVSLGHGSILEENVLLGSGTVIGSNCFITNS
VIGPGCHIGDNNVLDQTYLWQGVRAAGAQIHQSLLC
DNAEVKERVTLKPRSVLTSQVVVGPNITLPEGSVISLH
PPDAEEDEDDGEFSDDSGADQEKDKVKMKGYNPAEV
GAAGKGYLWKAAGMNMEEEEELQQNLWGLKINMEEE
SESESEQSMDSEEPDSRGGSPQMDDIKVFQNEVLGT
LQRGKEENISCDNLVLEINSLKYAYNVSLKEVMQVLSH
VVLEFPLQQMDSPLDSSRYCALLLPLLKAWSPVFRNYI
KRAADHLEALAAIEDFFLEHEALGISMAKVLMAFYQLEI
LAEETILSWFSQRDTTDKGQQLRKNQQLQRFIQWLKE
AEEESSEDD

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Entrez GeneID: 8893