

Datasheet

TBRG4 mouse monoclonal antibody (hybridoma)

Gene Alias: CPR2, FASTKD4, KIAA0948

Catalog Number: H00009238-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant TBRG4.

Immunogen: TBRG4 (NP_004740.2, 1 a.a. ~ 631 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MAAHLVKRCTCLLREARQAPAMAPVGRLRLAWVAH
KTLTSSATSPISHLPGSLMEPVEKERASTPYIEKQVDHL
IKKATRPEELLELLGGSHDLDSNQAAMVLIRLSHLLSEK
PEDKGLLIQDAHFHQLLCLLSQIASVWHGTL SKLLGS
LYALGIPKASKELQSVEQEVWRMRKLKYKHLAFLAE
SCATLSQEQHSQELLAELLTHLERRWTEIEDSHTLVTV
MMKVGHLSLEPLMNRLEDKCLELVEHFGPNELRKVLVM
LAAQSRRSVP LLRAISYHLVQKPFSLTKDVLLDVAYAY
GKLSFHQTQVSQRLATDLLSLMPSLTSGEVAHCAKSF
ALLKWL SLPLFEAFAQHVLNRAQDITLPHLC SVLLAFAR
LNFHPDQEDQFFSLVHEKLGSELPGLEPALQVDLVWA
LCVLQQAREAE LQAVLHPEFHIQFLGGKSQKDQNTFQ
KLLHINATALLEYPEYSGPLL PASAVAPGPSALDRKVTP
LQKELQETLKGLLG SADKGSLEVATQYGWVLD AEVLL
DSDGEFLPVRDFVAPHLAQPTGSQSPPPGSKRLAFLR
WEFPNFNSRSKDLLGRFVLARRHIVAAGFLIVDVPFYE
WLELKSEWQKGAYLKDKMRKAVAEELAK

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Entrez GeneID: 9238

Gene Symbol: TBRG4