

Datasheet

LRRC14 mouse monoclonal antibody (hybridoma)

Catalog Number: H00009684-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant LRRC14.

Immunogen: LRRC14 (NP_055480.1, 1 a.a. ~ 493 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MHTLVFLSTRQVLQCQPAACQALPLLPRELFPLLFKVA
FMDKKTIVLRELVTWPFPLLSFQQLQECACSRAL
LQERPSTESMQAVILGLTARLHTSEPGASTQPLCRKH
ALRVLDMTGLLDDGVEQDPGTMSMWDCTAAVARTCI
AQQQGGAEPGPAPIPEVRVVDLRVNRASYAFLREAL
RSSVGSPLRLCCRDRLRAEDLPMRNTVALLQLLDAGCL
RRVDLRFNNLGLRGLSVIIPHVARFQHLASRLHYVHG
DSRQPSVDGEDNFRYFLAQMGRTCLRELSMGSSLL
SGRLDQLLSTLQSPLELELAFCALLPEDLRFLARSPH
AAHLKKLDLSGNDLSGSQLAPFQGLLQASAATLLHLEL
TECQLADTQLLATLPILTQCASRLYLGLYGNPLSMAGL
KELLRDSVAQAELRTVVHPFPVDCYEGLPWPPASVL
LEASINEEKFARVEAELHQLLLASGRAHVLWTTDIYGR
LAADYFSL

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Entrez GeneID: 9684

Gene Symbol: LRRC14

Gene Alias: KIAA0014