

Datasheet

AGR2 monoclonal antibody (M04A), clone 7E9

Catalog Number: H00010551-M04A

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant AGR2.

Clone Name: 7E9

Immunogen: AGR2 (AAH15503.1, 1 a.a. ~ 175 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MEKIPVSAFLLLVALSYTLARDTTVKPGAKKDTKDSRP
KLPQTL SRGWG DQLIWTQTYEEALYKSKTSNKPLMIH
HLDECPHSQALKKVFAENKEIQKLAEQFVLLNLVYETT
DKHLSPDGQYVPRIMFVDP SLTVRADITGRYSNRLYAY
EPADTALLLDNMKKALKLLKTEL

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In ascites fluid

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10551

Gene Symbol: AGR2

Gene Alias: AG2, GOB-4, HAG-2, XAG-2