

Datasheet

HTATIP2 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00010553-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human HTATIP2 protein.

Immunogen: HTATIP2 (NP_006401.2, 1 a.a. ~ 242 a.a) full-length human protein.

Sequence:

MAETEALSKLREDFRMQNKSVFILGASGETGRVLLKEI
LEQGLFSKVTLIGRRKLTfDEEAYKNVNQEVVDFEKLD
DYASAFQGHdVGFCClGTTRGKAGAEgFVRVDRDYV
LKSAELAKAGGCKHFNLLSSKGADKSSNfLYLQVKGE
VEAKVEELKFDRYSVFRPGVLLCDRQESRPGEWLVRK
FFGSLPDSWARGHSVPVVTVVRAMLNNVVRPRDKQM
ELLENKAIHDLGKAHGSLKP

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10553

Gene Symbol: HTATIP2

Gene Alias: CC3, FLJ26963, SDR44U1, TIP30