



ZPACK12321007 TIN BOTTLE LIQ 5 L UN Y



SPECIFICATIONS

UN-certification :	3A1/Y/180
Material :	Tin
Neck Type :	DIN 45
Filling Opening :	Flexout
Height :	300,00 mm
Length :	125,00 mm
Wide :	164,00 mm

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343