

Datasheet

ENHO polyclonal antibody

Catalog Number: PAB15531

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of ENHO.

Immunogen: A synthetic peptide corresponding to amino acids 10-60 of human ENHO.

Host: Rabbit

Reactivity: Bovine, Human, Mouse, Rat

Applications: WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: 100% sequence identity with human and mouse protein, 88% sequence identity with bovine protein. This antibody is useful for Western blot analysis where a band can be seen at ~7.9 KDa.

Form: Liquid

Recommend Usage: Western Blot (2 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (30% glycerol, 0.1% sodium azide)

Storage Instruction: Store at 4°C for short term. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 375704

Gene Symbol: ENHO

Gene Alias: C9orf165, UNQ470