




## Data Sheet

Image 1	 <p>Immunohistochemistry with Human Skeletal Muscle lysate tissue at an antibody concentration of 5.0ug/ml using anti-TNNI1 antibody (ARP42117_P050)</p>
Image 2	<div><div><p><b>Human HepG2</b></p><p>70 kDa 60 kDa 48 kDa 36 kDa 21 kDa</p></div><div><p><b>WB Suggested Anti-TNNI1 Antibody Titration:</b> 0.2-1 ug/ml <b>Positive Control:</b> HepG2 cell lysate</p></div></div>
Image 3	<div><div><p><b>Human Kidney</b></p></div><div><p>Rabbit Anti-TNNI1 Antibody Catalog Number: ARP42117 Paraffin Embedded Tissue: Human Kidney Cellular Data: Epithelial cells of renal tubule Antibody Concentration: 4.0-8.0 ug/ml Magnification: 400X</p></div></div>

This product is for Research Use Only. Not for diagnostic, human, or veterinary use.  
Optimal conditions of its use should be determined by end users.